
OPTIKA[®]
S C I E N C E
I T A L Y

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

PRELIMINARY OPERATIONS

We congratulate you on choosing this machine which will allow you to carry out your work with the greatest satisfaction and in complete safety.

THIS BOOKLET IS AN INTEGRAL PART OF THE MACHINE AND MUST BE READ CAREFULLY BEFORE INSTALLING, STARTING AND USING THE MACHINE.

The booklet contains important information and instructions regarding the SAFETY OF USE and MAINTENANCE of the pump and therefore must be stored with care.

DESCRIPTION AND INTENDED USE

The AV series vacuum pumps are rotary vane pumps, lubricated by oil injection. The coupling is direct without the interposition of a joint in order to obtain maximum compactness. The suction fitting incorporates a filter (to prevent foreign bodies from being sucked in) and a check valve to prevent possible oil leaks into the drained system. To be sure of no oil rising in the system, it is possible, upon request, to mount a solenoid valve on the suction fitting directly connected to the electric motor.

These pumps are equipped with a ballast valve that allows the purification of the oil in the event that condensable vapors must or have been aspirated.

As regards the electrical part, the motor has been built according to the European safety standards and complies with the IEC 34 (EN 60034) standard and the EEC 73/23 (1973), EEC 93/68 (1993), EMC EEC 89 directives. / 336 and EEC 89/392. It also has a thermal protector with an intervention threshold of 130 ° C inside to avoid any damage to the winding in any eventuality (for example pump seizure, motor failure, power supply voltage too low / high, excessive overheating, etc ...).

The AV series pumps are intended exclusively to create vacuum in tanks, systems, chambers, ovens, etc ...

They can be conducted by hand or installed on equipment or frames; in the latter case, the fixing of the pump must be performed in a workmanlike manner and must not present any danger to things and / or people.

THE INTENDED USE OF THE MACHINE MUST BE SCRUPULOUSLY RESPECTED. ANY OTHER USE IS TO BE CONSIDERED IMPROPER.

THE MANUFACTURER CANNOT BE CONSIDERED RESPONSIBLE FOR DAMAGES CAUSED BY IMPROPER USE OF THE MACHINE OR INCORRECT INSTALLATIONS OF THE SAME.

THE MACHINE MUST NOT BE TAMPERED FOR ANY REASON.

OTHERWISE, THE MANUFACTURER DECLINES ALL LIABILITY FOR THE OPERATION AND SAFETY OF THE MACHINE.

PRELIMINARY OPERATIONS

Unpack the pump and make sure it is intact and undamaged. Otherwise, do not use the machine and contact your dealer. For transport reasons, THE PUMP IS BOXED WITHOUT LUBRICANT. A bottle of oil sufficient for the first charge is provided in the packaging. For this operation, follow the instructions contained in the "installation" chapter.

Make sure that the packaging elements (bags, boxes, polystyrene sides) are placed out of the reach of children. Make sure, before use, that the pump is provided with plates (one on the motor and one on the carter of the machine). Otherwise, do not use the machine and notify the dealer immediately. The plate placed on the motor lists all the data necessary for a correct connection to the electrical network. The one placed on the pump casing indicates the technical data of the pump.

N.B. The removal of these labels may cause damage and / or inconvenience, as well as voiding the warranty.

INSTALLATION

THE MACHINE MUST BE PLACED ON A SAFE PLANE, IN A HORIZONTAL POSITION.

Oil charge

Use the bottle of oil supplied with the pump to carry out the first charge. Unscrew the breather cap (1) or the oil separator filter (if the machine is equipped with it) and pour in the oil until the level indicated by the warning light (4) is reached, which must be halfway through it. Slightly higher quantities do not cause malfunctions. Run the pump with its mouth closed disconnected from the system for a couple of minutes and then check the level again; top up if necessary.

Suction connection

The pump must be connected to the system through the use of neoprene, vacuum para rubber or reinforced PVC pipes; other materials such as copper or rilsan are also allowed as long as suitable fittings are used. The minimum diameter of these pipes is 10 mm and you have to try to reduce their length to a minimum and assign them the most straight path possible. Lower section measurements and less direct paths would only (negatively) affect the performance of the pump. If rigid pipes are used, it is necessary to avoid that they have to force on the pump connection.

Steam exhaust

The pump can operate without drain fittings, however, in the event that large capacities have to be emptied or that emptying is continuous and frequent, we recommend replacing the standard breather plug with a specially made oil separator filter. It is also possible to convey these vapors through pipes provided that they have a minimum internal diameter of 20 mm and that they do not cause impediments to discharge into the atmosphere.

DO NOT OBSTRUCT OR PLUG THE PUMP VENT HOLE AS THIS COULD CHECK DANGEROUS OVERPRESSURE INSIDE THE CRANKCASE!

INSTALLATION

Connection to the electrical network

Always make sure that the mains voltage corresponds to that indicated on the identification plate on the motor.

Always make sure that the pump is connected to an electrical system with an efficient earth line.

Make sure that the socket is protected with a "circuit breaker" (differential magnetothermic switch with sensitivity less than 30 mA).

Insert the plug only after making sure that the machine switch is in the OFF (O) position.

If the pump is driven by a three-phase motor, check that the direction of rotation coincides with the direction of the arrow on the motor fan cover.

**IF THE CURRENT GOES OUT DURING USE,
SWITCH OFF THE MACHINE FOR SAFETY REASONS.**

STARTING UP

Each pump is scrupulously tested and is subjected to a preliminary running-in in the factory; however, it is good practice, at the first start-up or after a prolonged period of inactivity (2 ÷ 3 months), to observe the following rules:

Make sure that the lubricating oil has been loaded using the level sight glass (4) located on the side of the crankcase.

Put the pump into operation by letting it run with its mouth closed for a few minutes, plugging the intake (5) with a piece of rubber of adequate size. This operation has the purpose of ensuring the lubrication of the internal parts of the pump.

Check the oil level with the pump running.

For subsequent start-ups it will be sufficient to check the oil level. Should any irregularities be found, it is always advisable to stop the pump and ask for instructions from qualified personnel or contact our assistance service directly.

STOPPING THE PUMP

Stopping the machine requires no special attention except to put it in the air (have it sucked directly from the atmosphere); this is useful to facilitate the subsequent restart of the pump.

If you want to preserve the vacuum in the chamber, you need to insert a valve that isolates the system from the pump.

On request, a solenoid valve is available as an accessory which, being connected directly to the motor, automatically isolates everything upstream from the pump (even in the event of accidental power failure). For repetitive cycles at short intervals it is advisable not to stop the pump to avoid excessively straining the motor. Before stopping the pump (eg at the end of the work day) it is advisable to run it for a few minutes with the mouth closed by opening the ballast valve in order to purify the oil. Make sure that the machine is placed out of the reach of children and unauthorized personnel. While not reaching very high temperatures at full speed, check for small burns) or things. The machine is built in compliance with the safety criteria prescribed by the regulations in force. The use of electrical appliances, however, involves the observance of some elementary rules.

Do not come into contact with live electrical parts; inspection, maintenance and repair operations must be carried out by qualified personnel. In any case, disconnect the plug before proceeding with the aforementioned operations.

The electric cable of the motor is of the molded type with schuko plug and complies with the electrical and safety standards currently in force. If you want to modify it (e.g. replace the plug) make sure to proceed according to the regulations in force.

If the pump is installed on equipment or systems, make sure that the fastening is safe and does not give rise to any danger to the users. Furthermore, since it is an air-cooled machine, check that the motor cover is spaced at least 4 cm from any blind part and in a sufficiently ventilated place.

Try to avoid directly exposing the machine to atmospheric agents. For example, do not expose the pump directly to the sun's rays in summer to avoid overheating the engine excessively, or to avoid winter frosts in order not to have starting problems.

Considering that the pump uses vacuum to suck the oil from the crankcase and lubricate its internal parts, avoid working at pressures above 100 mbar for prolonged periods.


For safety purposes, use only original accessories and spare parts.

THE MANUFACTURER DECLINES ANY LIABILITY RELATING TO ANY DAMAGE THAT MAY RESULT FROM FAILURE TO COMPLY WITH THE INTENDED USE, INSTRUCTIONS AND WARNINGS SHOWN IN THIS BOOKLET.

ORDINARY MAINTENANCE

BALLAST VALVE (GAS BALLAST)


The AV-12 series pumps are equipped with a ballast valve to open in the event that condensable vapors need to be sucked. The suction of these vapors would pollute the pump lubricant, preventing it from reaching the ultimate vacuum and, in extreme cases, could compromise its pump operation. The opening of this valve, obtained by unscrewing the knob a few turns, has the function of decreasing the compression ratio inside the pump and preventing the condensation of these vapors which, precipitating, would pollute the oil. Once these vapors have been aspirated and discharged into the atmosphere, close the gas ballast to reach the maximum degree of vacuum.

OPTIONAL SEPARATOR FILTER

If large capacities have to be emptied or if you have to do continuous and frequent emptying, it is advisable to install an oil separator filter specially made by us. This accessory, supplied on request, replaces the standard breather plug (1) and has the function of avoiding the smokiness at the exhaust (micro-drops of oil carried by the air flow during passage in the pump chamber) which occurs during the first phases of emptying. The filter cartridge separates the air from the oil drops by coalescence collecting them in the filter holder. The only maintenance to be carried out is that of

frequently check that too much oil is not collected inside the filter body and eventually drain it using the appropriate cap. The life of the cartridge varies between 1,500 and 2,500 hours depending on the service to which it is subjected.

Correct observance of the rules contained in this paragraph is the basis for a long life of the pump and for achieving the declared performance.

Types of oil to be used.

A bottle of oil sufficient for the first charge is supplied with the pump. For subsequent changes, the oil to be used must be of a mineral type without additives with a viscosity of 46 cSt. With the DIP type pumps it is also possible to use synthetic oil based on POE polyol esters, in case it is necessary to treat HFC refrigerant gases (R134a, R404A, etc ...). Below is a table with some types of compatible mineral oil (highlighted is the first supply oil).

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AGIP	MAG 46
CASTROL	PERFECTO GG
ELF	MOVIXA 46
ESSO	NURAY 46
IP	HADRIA OIL 38
MOBIL	RUBEX 300
SHELL	CARNEA OIL 46

ORDINARY MAINTENANCE

Oil change

Periodically check the oil level and its degree of purity through the indicator light (5). The first oil charge must be replaced after about 150 hours of operation; the subsequent ones, if no pollution phenomena occur, after about 700 ÷ 900 hours and in any case at least every 6 ÷ 7 months. To replace the oil, proceed as follows: drain the oil contained in the pump casing using the appropriate cap, preferably when the machine is hot. To obtain a complete emptying, after draining, run the pump with its mouth open for about 30 seconds, then let the oil drain. Screw the drain plug back on and charge as described in the installation chapter.

EXTRAORDINARY MAINTENANCE

AS LUBRICATING OIL IS A POLLUTING PRODUCT, IT MUST NOT BE DISPOSED OF IN THE ENVIRONMENT, IT IS A SPECIAL WASTE AND HOW SUCH IT MUST BE DISPOSED OF ACCORDING TO THE LAWS IN FORCE!

DURING OIL UNLOADING AND LOADING OPERATIONS USE SUITABLE PROTECTIVE GLOVES.

It must be carried out when serious anomalies occur, such as overheating, significant vacuum drops, abnormal noise or pump blockage resulting from the aspiration of foreign bodies. In these cases, disassembly of the pump, careful washing and appropriate replacements or repairs are certainly required. This type of intervention must be entrusted to qualified personnel following the rules set out below and the exploded drawing attached to the instructions.

Referring to the identification numbers and the bill of materials, the necessary spare parts may be requested. In any case, as the user does not have suitable personnel, he can always contact our after-sales service for a perfect set-up of the pump.

Rules for disassembly

To disassemble the pump, proceed by following the instructions listed below after disconnecting the motor from the electrical system and draining the oil.

- Remove the pump casing by unscrewing the four fixing screws.
- Proceed to disassemble the front cover by unscrewing the three nuts.
- Remove the pump body.
- With the aid of an extractor remove the rotor from the shaft and collect the vanes.
- Disassemble the radial compensator, taking care not to deform it excessively, and the safety spring pin.
- Remove the rear cover and its gasket.

EXTRAORDINARY MAINTENANCE

If it is a single-stage pump, the operation is complete; if, on the other hand, it is a two-stage model, before removing the rear cover, proceed as follows:

- Remove the diaphragm.
- Detach the second stage body from the intermediate flange.
- Using an extractor, remove the rotator from the shaft and collect the vanes.
- Remove the radial compensator ring (always taking care not to deform it) and the safety pin.
- Remove the rear cover and its gasket.
- At this point the competent mechanic is able to assess the cause that led to the intervention and to proceed with any replacements and / or repairs. If, on the other hand, the cause of the imperfect operation should be attributable to a loss of the sealing ring placed in the pump flange, it will be necessary to disassemble the motor. To do this, you need to:
 - Remove the fan cover by unscrewing its four screws.
 - Remove the cooling fan from the motor shaft.
 - Unscrew the four motor fastening tie rods and remove the motor casing.
 - After removing the shaft from the flange, the gaskets can be replaced, taking care, when reassembling them, to position them in the same direction as the previous ones.
 - To reassemble the motor and pump, proceed in the opposite direction to what was done for disassembly, carefully cleaning and degreasing the support surfaces and all the disassembled parts. Pay particular attention to:
 - Do not damage the sealing ring by inserting the shaft into the flange.
 - Do not axially force the shaft bearings excessively during reassembly.
 - Make sure that the motor casing fits snugly against the coupling flange.
 - Take care to place the vanes in the rotor slots with the outer rounded surface.

THE SEAL BETWEEN THE PLANS OF THE PUMP BODIES AND THE LIDS IS OBTAINED USING LOCTITE 573 RESIN, THEREFORE YOU MUST OBTAIN THIS MATERIAL BEFORE PROCEEDING WITH THE INTERVENTION.

- THE PUMP CRANKCASE GASKET MUST BE REPLACED WHEN THE CRANKCASE IS DISASSEMBLED.
- DISASSEMBLY OF THE PUMP BY THE CUSTOMER AUTOMATICALLY VOIDS THE WARRANTY.


GENERAL WARRANTY CONDITIONS

The AV series vacuum pumps are covered by a full mechanical-electrical guarantee for a period of 24 (twenty-four) months from the date of purchase. Any damage caused by non-compliance with the use and maintenance rules mentioned in this instruction booklet cannot fall under the pump warranty.

Each machine to be repaired (whether under warranty or not) must be delivered ex works. factory without lubricating oil and possibly accompanied by a brief indication of the anomaly found.

ECOLOGICAL MEASURES

The following symbol is applied to the device:

This symbol indicates that the device can be a source of precious raw materials, and therefore must not be thrown away, but collected separately and sent for recovery and recycling. These raw materials, if dispersed in the environment, they can cause damage to the environment itself and to the health of citizens.

The device manufacturer, M.A.D. Scientific Equipment - Optika Microscopes, as well as the Public Administration, undertake to protect the environment by recovering, reusing and recycling the raw materials contained in the device. However, this recovery requires your contribution.

When, at the end of its use, the appliance has to be disposed of, do not try to open it, reuse parts in a different way from that described in this User Manual, and do not throw the appliance in the trash, but return it to the dealer from the which you have purchased, or to another M.A.D. od Optika Microscopes. The retailer will collect the device free of charge.

The recovery of the raw materials will therefore take place according to the European Directive 2002/96 / EC, its updates and all other Directives in force.

Attention: the law provides for penalties for those who illegally dispose of equipment bearing the above symbol. M.A.D. Scientific Equipment - Optika Microscopes disclaims any liability to anyone who disposes or modifies the equipment in a way that differs from what is indicated here.

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OPTIKA S.r.l - Copyright